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Table 1. MW measurement repeatability results

Sample ID	Reference MW, kg/mol	Repeat	MW by NMR, kg/mol	Average MW for 5 repeat NMR measurements, kg/mol	Standard deviation for 5 repeat NMR measurements, kg/mol	Relative standard deviation for 5 repeat NMR measurements, %
Sample 1	864	1	803	800	24	3
		2	786			
		3	771			
		4	844			
		5	797			
Sample 2	8980	1	8548	8543	179	2
		2	8660			
		3	8554			
		4	8215			
		5	8738			


Figure 1. Calibration for measuring MW in polyolefin samples. Correlation coefficient  $R^2=0.9924$ , standard deviation 0.0956, variance 0.0091.

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This publication is the copyright of Oxford Instruments plc and provides outline information only, which (unless agreed by the company in writing) may not be used, applied or reproduced for any purpose or form part of any order or contract or regarded as the representation relating to the products or services concerned. Oxford Instruments' policy is one of continued improvement. The company reserves the right to alter, without notice the specification, design or conditions of supply of any product or service. Oxford Instruments acknowledges all trademarks and registrations. © Oxford Instruments plc, 2019. All rights reserved. Part no: MR/217/0619

The Business of Science®